

CUSTOMER NOTICE

1. Introduction

The purpose of a Customer Notice is to provide timely communication regarding updates to Leica Biosystems product user manuals. This Customer Notice contains information about your Leica Biosystems product that will be included in a future release of one or more user manuals for your product. This information may include new instructional content, updated specifications, or additional information about the use or care of your product. Use this document as supplementary information until the revised user manuals are available.

This Customer Notice provides information and instructions for using the new features provided with the Aperio GT 450 Release 1.3.1.

This document contains the following sections:

1. Introduction
2. Aperio GT 450 Release 1.3.1 overview
3. Set the Pyramid Level Downscaling option in SAM

Products Affected

Aperio GT 450 Scanner, Aperio GT 450 Scanner Administration Manager Software, and Aperio GT 450 Controller Version 1.3.1.

User Manuals Affected

The information contained in this document will be incorporated into the following user manuals upon their next major scheduled release:

Document Title	Base Document Number
Aperio GT 450 IT Manager and Lab Administrator's Guide	MAN-0394

2. Aperio GT 450 Release 1.3.1 overview


The Aperio GT 450 release 1.3.1 includes the following feature updates:

2x pyramid level downscaling support

With Aperio GT 450 Release 1.3.1, you have the option of changing the pyramid level downscaling from the default of 4x to 2x. This option is set in the Aperio GT 450 SAM software.

When displaying Aperio GT 450 scanned images, some third-party viewing software solutions must interpolate the image data as the user navigates through the available magnification levels. This interpolation process can cause a decrease in viewing performance. Setting the pyramid level downscaling to 2x can increase viewing performance by enabling third-party viewing software solutions to display all available magnification levels with less interpolation of the image data.

The illustrations below show how 2x pyramid level downscaling provides the viewing software direct access to a greater range of magnification levels.



When you enable the 2x pyramid level downscaling option, the image size for 2x downsampled images is approximately 15% to 20% greater than the default 4x downsampled images, depending on the slide image.


For more information on enabling 2x pyramid level downscaling in SAM, see [Set the Pyramid Level Downscaling option in SAM \(on page 3\)](#).

3. Set the Pyramid Level Downscaling option in SAM


Only a user who is assigned the Lab Admin role can make configuration changes.

To adjust the new Pyramid Level Downscaling settings for a scanner, follow these steps:

1. Log in to the SAM software.
2. When the SAM Home screen appears, confirm that the Scanners icon in the banner is selected, and the page shows the list of scanners. Click the **Scanners** icon to display the list, if necessary.



3. Click the **System Information** icon for the appropriate scanner.



4. In the left pane of the System Information page, click **Settings** to open the Scanner Config page.

Scanners Users

Aperio GT 450 SAM Software | UDI: 00815477020471(8012) 1.3
LeicaAdmin

SS45073 GT450

System Information Event Logs Configuration ONLINE

Info Scanner Config

Scanner Statistics Camera Config

Settings Scanner Additional Config

Focus Algorithm Config

RT Camera Config

RT Focus Config

Tissue Finder Config

Motion Config

Autoloader Config

https://10.221.226.116/#

Scanner Config

MACROFOCUS START

12.588305

MACROFOCUS END

11.588305

MACROFOCUS RESOLUTION

0.000125

MACROFOCUS RAMPOUT

0.1

MACROFOCUS POS OFFSET

0

Export XML Import XML Edit

- On the Scanner Config page, click **Edit**, and then scroll to the bottom of the page.
- In the **PYRAMID LEVEL DOWNSCALING** field, select **2x** to enable the 2x Pyramid Level Downscaling option.

Scanners Users

Aperio GT 450 SAM Software | UDI: 00815477020471(8012) 1.3
LeicaAdmin

SS45073 GT450

System Information Event Logs Configuration ONLINE

Info Scanner Config

Scanner Statistics Camera Config

Settings Scanner Additional Config

Focus Algorithm Config

RT Camera Config

RT Focus Config

Tissue Finder Config

Motion Config

Autoloader Config

Debug Options

COMPRESSION NUM THREADS

8

COMPRESSION QUALITY

91

COMPRESSION TILE SIZE

256

ICC PROFILE PATH

Segue_v3.icc

IMAGE FILTER

3

GAMMA

2.2

MINIMUM SLIDE WIDTH

22

PYRAMID LEVEL DOWNSCALING

2x

2x

4x

Export XML Import XML Save Cancel

- Click **Save**.

Notices

CN-0004, Revision A | February 2024

This manual applies to Aperio GT 450 Controller, Aperio GT 450 Console, and Aperio GT 450 SAM versions 1.3.1 and later.
Original Instructions.


Copyright Notice

- Copyright © 2024 Leica Biosystems Imaging, Inc. All Rights Reserved. LEICA and the Leica logo are registered trademarks of Leica Microsystems IR GmbH. , GT, and GT 450 are trademarks of Leica Biosystems Imaging, Inc. in the USA and optionally in other countries. Other logos, products, and/or company names might be trademarks of their respective owners.
- This product is protected by registered patents. For a list of patents, contact Leica Biosystems.

Customer Resources

- For the latest information on Leica Biosystems products and services, please visit [LeicaBiosystems.com/Aperio](https://www.leicabiosystems.com/Aperio).

Contact Information – Leica Biosystems Imaging, Inc.

Headquarters	Customer Support
 <p>Leica Biosystems Imaging, Inc. 1360 Park Center Drive Vista, CA 92081 USA Tel: +1 (866) 478-4111 (toll free) Direct International Tel: +1 (760) 539-1100</p>	<p>Contact your local support representative with any query and service request.</p> <p>https://www.leicabiosystems.com/contact-us/</p>

For research use only. Not for use in diagnostic procedures.



UDI 815477020228, 815477020464, 815477020471, 815477020563

REF 23GT450, 23GT450SAM, 23GT450SAMSW, 23GT450ZSTACKSW