

ELEVATE YOUR DIAGNOSTIC EXPERIENCE

APERIO eSLIDE MANAGER WITH WEBVIEWER AND LIS CONNECTIVITY



For In Vitro Diagnostic Use. The clinical use claims described for the products in the information supplied have not been cleared or approved by the U.S. FDA or are not available in the United States.

Advancing Cancer Diagnostics
Improving Lives

Leica
BIO SYSTEMS

APERIO eSLIDE MANAGER WITH WEBVIEWER AND LIS CONNECTIVITY

The number of pathology service requests arriving each week is rising, while the number of properly trained professionals is barely keeping pace. The tools your diagnostic laboratory needs to stay ahead of this increasing workload must evolve beyond optical instruments invented 400 years ago. **Fast forward to your future.**

Begin your digital pathology journey with high-quality whole slide images, obtained from high-throughput Aperio scanners, and retain precision quality with Aperio eSlide Manager image management software with WebViewer interface. Designed to optimize digital pathology processes for pathologists, laboratory professionals and IT managers, this high-quality image management solution is integrated directly into your existing workflow with the Aperio LIS Connectivity module.

Aperio software elevates your diagnostic experience using the power of digital pathology – manifesting the potential of universal case access and borderless professional collaboration.




ENRICH YOUR IMAGE QUALITY

START WITH QUALITY – STAY WITH QUALITY
TO PROMOTE BETTER RESULTS

High quality digital scans deserve an image management system that detects, stores and delivers whole slide images and metadata exactly where and when needed for diagnosis:

- Synchronizations can be done with past and current cases at the same time, in the same window.
- The conference functionality within the Aperio eSlide Manager platform enables pathologists to collaborate and share their knowledge and medical opinions.




ENHANCE YOUR IMAGE MANAGEMENT

METADATA CAPTURE
SUPPORTS OPERATIONAL EFFICIENCY

The LIS worklist is your entry point into an Aperio curated case space:

- Aperio eSlide Manager works in the background, using barcoded metadata to auto-assemble Whole Slide Images into LIS-designated cases.
- WebViewer launches directly from the LIS based on the metadata in Aperio eSlide Manager.
- Advanced searching of slide metadata facilitates case management with retrieval of images from past cases, for comparison purposes.



EXTEND LIS-DRIVEN WORKFLOW

A **PRACTICAL DIGITAL PATHOLOGY WORKFLOW**
 REQUIRES REAL-TIME, VALUE-ADD LIS INTEGRATION

The LIS Connectivity module unites our Aperio eSlide Manager workflow with your specific LIS-based administrative process to improve the pathologist's diagnostic experience.

- Once Aperio eSlide Manager has been integrated with your LIS, the pathologist is able to operate from an LIS worklist and launch WebViewer directly to consolidate their "case space."
- Both systems are kept in sync with industry standard HL7 Communication protocols.



* The LIS Connectivity module is available for Aperio eSlide Manager, but it is a separate product and additional costs are associated with it.

For In Vitro Diagnostic Use. The clinical use claims described for the products in the information supplied have not been cleared or approved by the U.S. FDA or are not available in the United States.

APERIO eSLIDE MANAGER WITH WEBVIEWER AND LIS CONNECTIVITY


YOUR FUTURE IN FOCUS

We have a clear vision of future pathology workstreams which will be aided by a wealth of artificial intelligence applications, but we know our success is about helping laboratories execute brilliantly now.

A COMPREHENSIVE DATA MANAGEMENT SYSTEM THAT YOU CONTROL

- Easily reconfigure the default data fields to capture valuable metadata unique to your organization and customize page layouts
- Customized configurations of metadata and streamlined workflow that you create is retained even through software updates
- Tag any data field that you might use later as a filter and your savable search criteria becomes as unique as you need it to be
- Ensure that the right images get to the correct people with a robust user permission structure

Consolidate Your Case Space





ESSENTIAL IMAGE-ENHANCED WORKFLOW FOR PATHOLOGISTS

The Aperio eSlide Manager WebViewer serves the pathologist as both a browser-based digital whole slide image viewer and as a portal into the Aperio eSlide Manager, the epicenter of the Aperio digital pathology solution. This dynamic, diagnostic tool securely connects to a gateway Data Server that controls access to metadata and scanned images located on a scalable Image Server infrastructure.

Advantages for Pathologists:

CONFIDENCE

- Compare current case with past patient cases to provide needed context, with the ability to synchronize image magnification and rotation of multiple images within the same window.
- Corroborate your diagnosis with colleagues and consulting sub-specialists using advanced communication features.
- Streamline slide navigation for professional quality presentations using handy Pin Drops.

EFFICIENCY AND FLEXIBILITY

- WebViewer launches directly from the LIS worklist based on the metadata in Aperio eSlide Manager
- Variable magnification (from image fit to highest scan resolution) and 360° rotation
- Annotations delineate areas for a special stain or a harvest zone for ancillary studies



ONE VENDOR, TOTAL SOLUTION FOR LABORATORY PROFESSIONALS

The advantages of working with one end-to-end digital pathology partner increase as the number of slides processed per day ramps up. The Aperio consulting and implementation services can help you navigate even the most complex deployments involving multiple system integrations and components, ensuring you have a single point of contact and resolution.

Advantages for the Laboratory Professional:

QUALITY

- Maintain your lab as a center of excellence in digital pathology

EFFICIENCY

- Streamlined workflow keeps histotechnicians and pathologists motivated and productive
- Keep track of metadata created in histology procedures such as grossing and embedding

REPUTATION

- Sustain an exemplary reputation for your lab in the areas of high-quality diagnoses and turnaround time

EASY TO MANAGE

- One vendor helps you manage the entire digital pathology workflow – from specimen capture to final report
- LIS integration keeps lab staff and pathologists happy and efficient, gives you peace of mind, and reduces redundant data entry



SECURE, SCALABLE ARCHITECTURE FOR IT MANAGERS

Aperio eSlide Manager offers a fast, secure, and flexible IT architecture. Scanner hardware can be near histotechnicians and image servers near case reviewers – all connected to a flexible image repository supporting a center of pathology excellence and teaching. Aperio eSlide Manager is also a proven robust data management solution that can handle hundreds of thousands of multi-gigabit image files and their metadata, as well as tracking and responding to a heavy stream of event messages.

Advantages for the IT manager:

SCALABLE

- Hub and spoke network architecture
- Modular application design
- Supports 150 concurrent users simultaneously accessing the database (Aperio eSlide Manager Performance Pack)

SECURE

- Single gateway for users across multiple locations
- Enhanced cybersecurity: PHP v7, updated Apache, OpenSSL for Data Encryption between browser and system web pages, and complex passwords
- In-app user authentication and authorization ensure that the right people have access to the right functionality and data

Server Requirements

Component	Server Minimum Requirements	Server Recommended Requirements	Performance Pack Server Requirements	Specific Reference Hardware*
CPU speed	Two dual-core or one Quad-core Intel 2.4GHz or faster	8-core processors 2.0GHz or faster recommended	CPU 1 Recent Generation 12-Core 2.3 GHz CPU	Intel Xeon Gold 5118 2.3G, 12C/24T, 10.4GT
Hard disk space	1TB DAS, NAS, or SAN Storage (recommended) RAID for redundancy	5TB DAS, NAS, or SAN Storage (recommended), RAID for redundancy	RAM 16 GB	2x 8GB RDIMM, 2666MTs
Available C: or OS drive space	100GB	100GB	Hard Drive 2x 240GB Mirrored SSD Hard Drives for OS and applications	Dell BOSS controller card + with 2M.2 Sticks 240G
Memory	4GB RAM	16GB RAM	Network Card 10GB	Broadcom 57416 2-Port 10Gb Base-T + 5720 2-Port 1Gb
Network card	1GB network card	Dual 1GB network cards	OS Windows Server 2016	OS Windows Server 2016
Power supply	Redundant power supplies (recommended)	Redundant power supplies		
Media input	Not required	DVD (optional)		
Operating system	WIndows Server 2016 minimum	Windows Server 2016		

Monitor Specifications

Feature	Specification
Display technology	LED
Resolution	1920 x 1200 pixel minimum
Screen Size (measured diagonally)	609.7 mm / 24.0 inches minimum
Contrast Ratio	1000:1 minimum
Color Gamut	99% sRGB minimum
Display Mode	sRGB
Brightness (also called luminance)	180 cd/m ² minimum
Display Colors	16.7M colors minimum
Connectivity	Digital interface (such as HDMI, DisplayPort, or USB-C)
Response Time	15ms maximum
Viewing Angle	178 degrees minimum
Non-reflective coating (may be listed as Glare Screen)	Yes

Viewing Workstation

Feature	Specification
CPU speed	4 Core @ 3.6 GHz minimum
Hard disk space	80 GB free disk space minimum
Memory	8GB or more recommended
Network card	1Gigabit network card or faster
Video card	2GB, 64-bit, 32GB/sec, 256 Cores minimum
Browser	Firefox, Chrome
Operating System	Windows 10 Enterprise 64 bit

* Recommended when expanding to more than 25 concurrent users

** Requirements are optional and only necessary if you need to support concurrent users greater than 25 (and up to 150)

For In Vitro Diagnostic Use. The clinical use claims described for the products in the information supplied have not been cleared or approved by the U.S. FDA or are not available in the United States.

Leica Biosystems is an international company with a strong network of worldwide customer services. For detailed contact information on your nearest sales office or distributor please visit our website:

LeicaBiosystems.com

Leica Biosystems is a global leader in workflow solutions and automation. As the only company to own the workflow from biopsy to diagnosis, we are uniquely positioned to break down the barriers between each of these steps. Our mission of "Advancing Cancer Diagnostics, Improving Lives" is at the heart of our corporate culture. Our easy-to-use and consistently reliable offerings help improve workflow efficiency and diagnostic confidence. The company is represented in over 100 countries. It has manufacturing facilities in 9 countries, sales and service organizations in 19 countries, and an international network of dealers. The company is headquartered in Nussloch, Germany. Visit LeicaBiosystems.com for more information.