Advancing Cancer Diagnostics Improving Lives



# Fast forward to next generation imaging

## Introducing the **Aperio GT 450**

Scale up Digital Pathology operations by increasing throughput in high volume pathology labs.

For research use only. Not for use in diagnostic procedures.

\*International patent families including patent applications and granted patents.

More than 30 Patents SMARTER\*

### Be on time Fast scans Confidence No shipping continuous loading Be more productive Next generation imaging for Technicians Auto QC Fast-paced First-time scan success Cross-training Easy to use

#### Designed to deliver rapid results with confidence for histotechnicians.



Assign priority cases.

(~

M

Ŭ

Automated image quality check and slide calibration during each scan.

- No-touch continuous rack loading during scanning.
- 450 slide capacity packed into small size.
- Excellent image quality using Leica optics.







32 sec scan speed at 40x\*.

81 slides per hour at 40x\*.

Scan 120,000 slides per year operating the Aperio GT 450 for around 5 hours/day\*.

\* 15mm x 15mm scan area @ 40x; assumes 5 hr. scan time per day for 300 work days per year



Aperio GT 450 open face design allows you to continuously load racks while scanning is in progress to maximize productivity.

Featuring automatic image quality control, so you have image quality confidence before your scans reach the Pathologist.

Stats page allows you to view slides scanned/day/ month/year, average scan area and average scan speed.

For research use only. Not for use in diagnostic procedures.

Login controls Easy setup Security Store and forward DICOM or SVS compatible Next generation imaging for IT Professionals Flexibility Large storage Encryption

Designed to deliver fast, secure, and scalable IT architecture for IT Professionals.

"The IT architecture for the Aperio GT 450, including the SAM central hub, is a gigantic leap forward to scale and support digital pathology IT operations."

- Chris Khacherian, UCLA Program Manager, Service Delivery, DGIT | Digital Technology

A secure, flexible and fast IT solution to scale up your Digital Pathology operations.





Cyber security including encryption and login controls.

HL7 and PACS ready, for **integration** into LIS, LIMS or PACS.



Securely scan slides at a **spoke** location and store images at a **Hub** location for all sites to log in and view images.



Store and forward: store images temporarily then transfer during low IT traffic times.



Set up user-specific **log in** for each Aperio GT 450.





No more workstations needed. SAM server and software allow IT managers to set up and monitor GT 450 scanners via the network, including dashboard status views, PIN controls and log out times.

For research use only. Not for use in diagnostic procedures.



### Faster turnaround time High resolution Image quality Limited rescans Image quality Integrated work-flow Craftsmanship state-of-the-art Next generation imaging for Pathologists Peace of mind

940631810

GERMANY

EICA

Precision Ease of use DICOM

Designed to rapidly deliver images with excellent image quality for Pathologists.

#### "I am really impressed with the sharpness and contrast of the images. It is definitely next-generation image quality."

- W. Dean Wallace M.D., Professor, Pathology and Laboratory Medicine, David Geffen School of Medicine at UCLA

> Aperio GT 450 utilizes a high performance 40x objective manufactured by Leica, the company that has engineered world class optics since 1841.



Aperio GT 450 illuminates tissue from both sides of the slide resulting in a 99.5% accurate tissue finder that recognizes faintly-stained tissue while excluding:

- > Pen marks
- > Stain residue
- > Scratches
- > Dust



Aperio GT 450 applications:

- > Digital archiving
- > Drug development studies
- > Cancer research
- > Toxicology research
- > Basic cell biology research
- > Immunology research

"Good cellular details at high power and good focus throughout the tissue." – Sally Agersborg, MD, PhD, Medical Director, NeoGenomics Laboratories, Inc.





From the lab to the IT room the Aperio GT 450 was designed to scale up digital pathology operations for Histotechnicians, IT Professionals and Pathologists



"The Aperio GT 450 has demonstrated in a workflow study at our laboratories substantial increases in case throughput while dramatically reducing hands on time for our technicians, which frees up additional FTE capacity that we can operationally utilize to respond to increasing annual slide volumes."

- Julie Broccardo, Director of Anatomic Pathology Operations, NeoGenomics Laboratories, Inc.

Aperio GT 450 Summary of Specifications	
Maximum Slide Capacity	450
Rack capacity	15 racks of 30 slides, 15 racks of 20 slides
Racks compatible with Leica HistoCore SPECTRA Workstation	Yes
Continuous rack loading while scanning	Yes
Priority rack scanning	Yes
Magnification	40x
Resolution	0.26 µm/pixel at 40x
Scan speed (15mm x 15mm area at 40x)	32 sec
Sustained throughput (15mm x 15mm area at 40x)	81 slides/hr
Automated Image Quality Check	Yes
Slide sizes accepted	1x3
Field of View (FOV)	1mm
Bar code engines	1D and 2D symbology
Dimension	20.8" (52.83 cm) W x 23" (58.43 cm) L x 19.5" (49.53 cm) H
Weight	140 lbs (63.5 kg)

Part Number	Description
23GT450	Aperio GT 450 Scanner
23GT450SAM	Aperio GT450 Scanner Admin Manager (SAM) Server
23GT450SAMSW	Aperio GT450 Scanner Admin Manager (SAM) Applications Software
23Dell740XD21TB	Dell 740 XD image storage DSR 21 Terabytes
23Dell740XD42TB	Dell 740 XD image storage DSR 42 Terabytes
23Dell740XD85TB	Dell 740 XD image storage DSR 85 Terabytes
23Dell740XD138TB	Dell 740 XD image storage DSR 138 Terabytes
23VS100	Aperio viewing station (includes calibrated monitor, workstation, Aperio ImageScope viewer and optimized ICC file)

#### For research use only. Not for use in diagnostic procedures.

Copyright © 2019 Leica Biosystems Imaging, Inc. All Rights Reserved. LEICA and the Leica logo are registered trademarks of Leica Microsystems IR GmbH. Aperio is a registered trademark of Leica Biosystems Imaging, Inc. in the USA and optionally in other countries. GT and GT 450 are trademarks of Leica Biosystems Imaging, Inc. in the USA and optionally in other countries. Other logos, product and/or company names might be trademarks of their respective owners.

\*Parts shown are each sold separately.

### \_earn more at: LeicaBiosystems.com/AperioGT450

190560 Rev A • IMC-1198-REV-A 08/2019