


Advancing Cancer Diagnostics
Improving Lives



All-In-One Remote Live View & Desktop Scanner



APERIO LV1 IVD


The clinical use claims described for the products in the information supplied have not been cleared or approved by the U.S. FDA or are not available in the United States.

Remote live view and desktop whole slide scanning in one compact package

Aperio LVI IVD is the compact, economic solution for your diagnostic needs. Remote slides are made available for live view in under 15 seconds, invaluable for time-dependent applications such as intraoperative consultation. It doubles as a small-footprint digital pathology scanner, enabling rapid scanning of frozen section or permanent H&E slides, for diagnostic review and archiving. Aperio LVI IVD delivers a cost-effective solution for on screen diagnosis, helping you to get results faster for your patients.

VIEW SLIDES REMOTELY IN

15



SECONDS OR LESS*

*Requires at least a 50 megabit per second (mbps) internet connection using Team Viewer software (v10) and Aperio LVI IVD console software (v4.0) with one slide loaded in the slide tray. 15 second clock starts once the Team Viewer password is entered and "Log On" is selected.

Key Applications



Diagnosis: Review and interpret H&E slides on-screen from your laboratory or office.




Remote Access: Share and control slides live in real time.



Intraoperative Consultation: Review frozen section slides with the surgical site.

Up to 16 simultaneous live views to compare same or different slides

Change magnification



Remotely navigate in x, y or z planes

Switch between or add slides

Clone views
Take snapshots


Add annotations or notes

Real Time Live View

- » Remote access to slides in 15 seconds
- » Ideal for intraoperative consultation & other time-sensitive uses
- » Provide answers for your patients fast
- » Accurate fine control – just like a microscope
- » View slides at up to 63x* magnification
- » Switch easily to scanner mode for archiving

Whole Slide Scanning

- » Compact & economic desktop scanner
- » 4 slide capacity with barcode support
- » Scan & view at 2.5x, 10x, 20x or 40x*
- » Integrated with Aperio eSlide Manager
- » Easily scan up to 5 z-stack layers
- » Intuitive design and 1-touch scanning




*63x digital zoom for viewing only



Consultation: Rapidly review without traveling to the




Archiving: Scan, capture, and archive a permanent copy of your slides in minutes.



Education: Facilitate interactive teaching with live view and scanned slides.

The clinical use claims described for the products in the information supplied have not been cleared or approved by the U.S. FDA or are not available in the United States.

APERIO LV1 IVD



LEICA BIOSYSTEMS

Leica Biosystems (LeicaBiosystems.com) is a global leader in workflow solutions and automation, integrating each step in the workflow. As the only company to own the workflow from biopsy to diagnosis, we are uniquely positioned to break down the barriers between each of these steps. Our mission of "Advancing Cancer Diagnostics, Improving Lives" is at the heart of our corporate culture. Our easy-to-use and consistently reliable offerings help improve workflow efficiency and diagnostic confidence. The company is represented in over 100 countries and is headquartered in Nussloch, Germany.

Leica Biosystems – an international company with a strong network of world-wide customer services. For detailed contact information on your nearest sales office or distributor please visit our website: LeicaBiosystems.com.

VISION 24

- IMPROVED QUALITY
- INTEGRATED SOLUTIONS
- OPTIMIZED EFFICIENCY

LeicaBiosystems.com/AperioLV1IVD

Copyright © 2018 Leica Biosystems Imaging, Inc. All Rights Reserved. LEICA and the Leica logo are registered trademarks of Leica Microsystems IR GmbH. Aperio and EC Plan-Neofluar are trademarks of the Leica Biosystems group of companies in the USA and optionally in other countries. Other logos, product and/or company names might be trademarks of their respective owners.

18756 Rev A IMC-1194-REV-A 07/2018