


HISTOCORE WATER BATH

HISTOCORE SLIDE DRYER

MAKE YOUR SLIDES FLOATING AND DRYING EASIER



HistoCore Water Bath:

- Detachable water tray
- Black water tray and LED illumination help to easier identify floating sections
- Special treatment improves tray durability
- The OLED touch panel with an intuitive user interface controls both water bath and slide dryer



HistoCore Slide Dryer:

- Space-saving design, small footprint
- Up to 30 slides capacity
- Black rack provides good visibility and contrast
- A 45° angled and heated rack facilitates slide pick up and placement
- No need for an additional power socket, simply connect it to the HistoCore Water bath

HISTOCORE
WATER BATHHISTOCORE
SLIDE DRYER

FLEXIBLE LAYOUT OPTIONS ADAPT TO YOUR WORKFLOW NEEDS



Model Name	Order Number	Dimension (L x H x W)	Heating Temperature
HistoCore Water Bath	140607020C0	280x280x105 mm	Up to 60 °C
HistoCore Slide Dryer	140607010C0	200x280x98 mm	Up to 75 °C

HistoCore Slide Dryer as accessory

Advancing Cancer Diagnostic
Improving Lives

Leica
BIO SYSTEMS