
APERIO CS2

HIGHLY RELIABLE, DESKTOP DIGITAL PATHOLOGY SCANNER



For Research Use Only. Not For Use In Diagnostic Procedures.



QUICKLY AND RELIABLY SCAN, DIGITIZE AND SHARE WHOLE SLIDE IMAGES

Produce high-quality digital slides right from your desktop with the ultra-compact Aperio CS2 digital pathology scanner. With a five-slide capacity and high resolution magnification capabilities, including 20x and 40x the Aperio CS2 is a highly reliable workhorse for medium-volume sites. The intuitive interface and easy-to-use design consistently provides rapid creation of high resolution digital slides, with >98% first time scan success rate.

HIGH QUALITY IMAGES

The Aperio CS2 provides rapid and seamless creation of whole slide, digital slides yielding high scan quality and resolution – occupying only 1.6 square feet of bench-top space.

SHORT TIME-TO-VIEW

With a fast 1GB/second network interface, and real-time image compression, the Aperio CS2 scans and moves images in record time for remote review. Digital slides are ready for review in the Aperio eSlide tray in 90 seconds.¹

SHARE ANYTIME, ANYWHERE

The ideal tool for creating digital slides that can be rapidly shared with colleagues, whether they are down the hall or across the globe. The Aperio CS2 can help to improve turnaround time on slide review, via remote access and assessment in research applications.



Easily load and scan up to five 1x3 slides



Slides are on the Aperio eSlide tray in <90 seconds



Review difficult cases with peers or collaborate globally

APERIO CS2 SCANNER

Reliable, easy-to-use, and compact; spend less time learning how to use the Aperio CS2 and more time improving productivity.


- Easily load 5-slide tray and scan slides automatically for Brightfield applications
- Complete an eSlide capture in minutes, then view anytime, anywhere without waiting for post-scan image processing
- Receive better first time success rates with fewer rescans.
- · Benefit from the advanced small footprint by creating high quality eSlides within limited space
- Reduce turnaround time for access to critical results



 Slide tray eliminates contact with delicate glass slides and remains horizontal at all times

Key Pathology Applications

- Image analysis
- Cancer Research
- Tissue Microarrays
- Cell Biology
- Local and global sharing



LeicaBiosystems.com



APERIO CS2



Letea Blosystems is an international company with a strong network of worldwide customer services. For detailed contact information on your nearest sales office or distributor please visit our website: LeicaBiosystems.com

LeicaBiosystems.com/Aperio

For Research Use Only. Not For Use In Diagnostic Procedures.

Advancing Cancer Diagnostic Improving Lives



Leica Biosystems is a global leader in workflow solutions and automation. As the only company to own the workflow from biopsy to diagnosis, we are uniquely positioned to break down the barriers between each of these steps. Our mission of "Advancing Cancer Diagnostics, Improving Lives" is at the heart of our corporate culture. Our easy-to-use and consistently reliable offerings help improve workflow efficiency and diagnostic confidence. The company is represented in over 100 countries. It has manufacturing facilities in 9 countries, sales and service organizations in 19 countries, and an international network of dealers. The company is headquartered in Nussloch, Germany. Visit LeicaBiosystems.com for more information.